SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

	Gebrauchsanweisung	2
	Instructions for use	7
	Manuel d'utilisation	12
	Istruzioni per l'uso	17
	Instrucciones de uso	22
	Instrukcje użytkowania	27
	Instruktioner för användning	32
	Instruktioner til brug	37
	Bruksanvisning	42
	Käyttöohjeet	47
	Gebruiksaanwijzing	52
	Οδηγίες χρήσης	57
	Upute za uporabu	62
	Instrucțiuni de utilizare	67
	Naudojimo instrukcijos	72
	Návod k použití	77
	Návod na použitie	82
	Használati utasítás	87

SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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Produkt/ Produkttyp

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Instructions for use

Product/ Product type

Co/ Cr blanks for the fabrication of customised dental restorations

Material type

Cobalt/chrome alloy

Class IIa medical device

Co/ Cr alloy type 4 according to DIN EN ISO 22674

The product is free from nickel, cadmium, beryllium and lead in accordance with DIN EN ISO 22674.

Users

Professional users who produce customised dental restorations

Target group

Patients of all ages and genders who require customised dental restorations

Indication / Intended use

The product is a medical device intended for the manufacture of dental prostheses for temporary or long-term use, which are partially introduced into the human body by means of a clinical procedure, but are inserted into the teeth or attached to the tooth structure with suitable fixing material and are therefore classified as Class IIa.

Indication:

- anatomically reduced crown and bridge frameworks in the anterior and posterior region
- Bridges with up to 14 units or bridges with small cross-sections
- Primary and secondary telescopic crowns
- clamps, bars, fasteners
- Fully anatomical crowns and bridges in the posterior region
- Restorations with thin cross-sections that are exposed to very high loads
- Free-end bridges with a maximum of one pontic
- implant-supported superstructures
- Removable partial dentures

The product is a type 4 Co/Cr alloy in accordance with DIN EN ISO 22674 and therefore has no indication restrictions compared to milled or cast type 4 Co/Cr alloys.

Contraindication

Do not use in case of proven hypersensitivity to one or more metals contained in the alloy

Veneering ceramics/ glaze

Co/Cr veneering ceramic

SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Material properties/ Technical data

Composition:		Technical data (after final sintering):	
Co	59,5 - 69%	0.2% proof stress	480 MPa
Cr	26,5 - 30%	Percentage elongation at break	22%
Mo	4,5 - 7%	Percentage breakage constriction	16%
Mn	0 - 1%	E-modulus	178 GPa
Fe	0 - 1%	Tensile strength	864 MPa
Si	0 - 1%	Density	7.59 g/cm ³
C	0 - 0,35%	Corrosion resistance	< 200 µg/cm ²
other	< 1%	Tarnish resistance	Yes
		Vickers hardness	224 HV1
		Coefficient of thermal expansion	14.26 x 10 ⁻⁶ /K

Description of the

The product consists of blanks made of a biocompatible cobalt/chromium alloy for the manufacture of dental prostheses.

The production method in sintered metals avoids the disadvantages of cast blanks. These can exhibit inhomogeneities, crystalline structures and different hardness ranges.

After sintering, the product is characterised by a homogeneous structure. It is easy to mill and the wear on milling cutters is low. The product can be milled both dry and wet.

Processing

Scaffolding production

The desired restoration is milled from the selected blank. The blanks can be processed with all commercially available dental CAD/CAM milling machines.

The resulting milling dust is very fine. Before milling, please check whether your milling machine is suitable for processing this material. Please observe the safety instructions.

As the material has a firm consistency but can be milled with extremely stable edges, milling machines for zirconium or non-ferrous metals can be used with the appropriate milling strategy. The sintering shrinkage of the material is indicated as a factor on the blank and must be set in the parameters of the milling machine.

Milling

For restorations of up to 5 units, a sintering support is only required for highly curved restorations. A sintering support must be created for restorations with 6 or more pontics. Connect each, but at least every second pontic with the sintering support. The thickness of the sintering support should be 1.5 mm, the connector thickness 1.4 mm. Terminal pontics should always be provided with a connector. (See illustrations.)

SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Minimum thicknesses:

Wall thickness of individual caps	0.4 mm
Edge thickness of individual caps	0.2 mm
Wall thickness of bridges	0.5 mm
Edge thickness bridges	0.2 mm

Bridges in the posterior region:

Connector cross-section	9 mm ²
Pendant at tooth position	Maximum 1
Connector cross-section to trailer	12 mm ²

Bridges in the anterior region:

Connector cross-section	6 mm ²
-------------------------	-------------------

Disconnect

Once the milling process is complete, the objects are cut out of the blank.

Due to the connector thickness of 1.4 mm, it may happen that a last connector remains when the restoration is separated, which, due to the considerably thinner marginal or wall thicknesses compared to the 1.4 mm of the connector, leads to a fracture of the coping if the restoration is only held by one connector. We therefore recommend that you first cut all connectors around the restoration by 50% to create a predetermined breaking point. In this way, you can ensure that the cut connector breaks first and not the thinner-walled coping of the restoration.

Sintering

Finally, any remaining milling dust is carefully brushed off the objects and placed in the sintering dish. They must not touch each other or the sintering dish, as this would cause the objects to fuse during the sintering phase. For optimum sintering results, the objects should be encased up to just below the crown margins. There should be no sintering beads in the objects or in narrow interdental spaces. Pontics, bar parts, etc. must not be covered, but should protrude at the top to prevent heat build-up. Wet objects must be pre-dried e.g. 80°C/30min.

The heating, sintering and cooling process is usually fully automatic. Please refer to the operating instructions for the respective furnace used. After the oven has cooled down to around 50°C, the sintered parts are removed.

SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

- ▶ Heating speed 15°C/ min
- ▶ Final temperature approx. 1280°C
- ▶ Holding time at final temperature 60 min
- ▶ Cooling uncontrolled and from 800°C with compressed air in a closed oven

Please note that when the product is sintered, the restoration acquires a grey to brownish oxide colour. The restoration can be blasted with a blasting medium (recommended grain size 110 µm) at a blasting pressure of 2-3 bar and then further processed.

Before sintering, please check that the argon bottle is sufficiently filled, that the lines are free of leaks and that the sintering bell is clean and undamaged. ZrO₂ sintering beads turn grey during sintering cycles under argon. This effect is normal and desirable.

In the event of a faulty sintering cycle without or with too little argon, grey sintering beads turn white again. The restoration oxidises completely and significant scaling occurs. The resulting oxide is green to blue in colour. The restoration no longer fits properly and must be remade. Re-sintering is not possible as the oxide is also located within the structure. After a sintering cycle without argon, at least one empty cycle with argon must be carried out. After this, the sintering beads must be grey again. If this is not the case and the beads remain white, check the fill level of the argon bottle and whether there may be a mix-up of the sintering beads (AL₂O₃ instead of ZrO₂). Check all sealing surfaces of the sintering dish, sintering bell and sintering plate for contamination. If the problem persists, please contact the manufacturer of your furnace.

Firing the ceramic

In principle, all commercially available ceramic bonding materials can be used. Please observe the processing instructions of the respective ceramic manufacturer and the CTE range specified therein for compatibility.

The minimum thickness of the finished copings should not be less than 0.3 mm. It is recommended to blast and clean (steam off) the frameworks with at least 110 µm aluminium oxide at 3-4 bar; the oxide firing is not mandatory, but can optionally be carried out for 5 minutes at 980 °C under vacuum (cleaning firing). The framework must be blasted with 110 µm aluminium oxide at 3-4 bar to thoroughly remove the oxide layer. Then clean the framework with a steam jet. When using a ceramic bonder, please follow the processing instructions of the respective manufacturer.

Soldering

We recommend a cobalt-based solder for soldering. Parts made from the product should not be soldered together with gold or palladium solder. The product is also very suitable for laser welding.

Safety instructions

Warning: Contains Cobalt (Co). The dust generated during processing of this product may cause cancer if inhaled, may impair fertility and may probably cause genetic defects. It is essential to observe the following safety measures.

Always wear respiratory protection (filter class FFP3), tight-fitting safety goggles, protective gloves and protective clothing and always switch on extraction equipment with filter class Hepa H. Do not inhale dust. Avoid contact with skin, mouth, eyes and clothing. Do not eat or drink while working. Keep away from food and drink. Wash hands after use. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke.

Storage

No special storage conditions. Store in the original packaging.

SRN: DE-MF-000021120

Basis UDI-DI: : ++EWPDCMSLQ

CopraSintec K

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Waste disposal

Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations. Do not dispose of with household waste. Do not allow to enter waterways, groundwater or sewage system.

Note

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the associated base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

Explanation of the labelling on the packaging

Confirmation: The product complies with the applicable European directives

Catalogue number/ item number

Batch designation

Unique product identification

Number of products in the packaging

Do not reuse

Observe the safety instructions in the operating instructions and instructions for use

www.wp-ifu.de/cocr-k

Path to the electronic instructions for use

Instructions for use in paper form available within 7 days

Medical device

Date of manufacture

Manufacturer

Swiss authorised representative

Rx only

Caution: Federal law restricts this product to sale by or on the order of a licensed physician or dentist.

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