rs VIVO

RS VIVO Tray

INSTRUCTIONS FOR USE

rapidshape

PACKAGING UNITS

RS VIVO Tray, 1000g [RS006198]

DESCRIPTION

RS VIVO Tray is a light-curing resin for the additive production of dental impression trays. The material is optimized for wavelengths of 385 and 405 nm in layer thicknesses of 150 µm and can be used with 3D printers of following manufacturers:

- Rapid Shape [D-Series, PRO-Series, ONE]

COMPOSITION

Preparation made of acrylates and phosphinate.

INDICATION

Additive manufacturing of dental impression trays.

CONTRAINDICATIONS

Patient and user: Product should not be used in patients or by users who are known to be allergic to any of the ingredients. Otherwise the patient group is not restricted.

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USE

BUILD PROCESS

- **Prepare data set:** The orientation and type of support structures can be found in the system supplier's instructions for use.
- **Establish variable process parameters:** The appropriate process parameters are selected according to the specifications of the printer system supplier. Ensure that the correct material-specific set of parameters is selected (product, layer thickness, exposure time).
- **Mixing materials:** Stir material to the bottom of the bottle.
First filling: Pour material into the material reservoir of the unit.
Refill: Refill the material as needed and mix in the reservoir with a silicone spatula. Allow the reservoir with material to rest in the unit for 10 min.
- **Start build process.**

CLEANING AND DRYING

- **Clean part:** Separate the parts (with support structures) from the build platform with a spatula. The component is cleaned under flowing isopropanol for **max. 5 min.** If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol. Cleaning with RS Wash / PROWASH (Rapidshape) devices is conducted by selecting the corresponding material parameter set (120s pre-cleaning / 65s cleaning / 90s drying / Nr. Of uses: 10).
- **Dry part:** Blow component with compressed air and leave to dry for 30 min at room temperature or in a drying cabinet adjusted to a temperature of approx. 37 °C.

POST-CURING

- For post-curing, the components have to be cleaned and cooled to room temperature.

Used Printer	Curing Device	Post-curing
Rapidshape	Otoflash	2000 flashes each for top and bottom / Σ 4000 flashes ¹
	RS Cure / PROCURE	Parameter set RS VIVO Tray ²

¹ Recommendation: Curing under exclusion of oxygen prevents the formation of the inhibition layer.

² exposure time: 6 min / vacuum 50 mbar / without nitrogen / power 100 % / wavelength 375+415 nm

Ensure that the unit used is serviced and calibrated in accordance with the manufacturer's instructions in order to preserve its mechanical properties, color and biocompatibility.

- Remove support structure if necessary.
- For post-curing without oxygen connection, clean under flowing conditions with isopropanol for max. 2 min. If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol.
- Blow component with compressed air and leave to dry.
- **POST-PROCESSING COMPLETED**

FURTHER PROCESSING INSTRUCTIONS

Parts are subsequently processed and finalized using standard dental techniques.

WARNINGS

Any deviation from the processing process described here can jeopardize the biocompatibility, functionality and safety of the final product. Contains skin sensitizing substances (acrylates, phosphinate). Please always observe the respective valid safety data sheets. Applies to the uncured material: May cause an allergic skin reaction. Avoid contact with unpolymerized material. Wear protective gloves / protective clothing / eye protection. If on skin: Wash with plenty of soap and water. Avoid breathing dust / vapours / spray. Switch on extraction system. If skin irritation or rash occurs: Get medical advice / attention. The unpolymerized material is harmful to the environment. Avoid release to the environment. Collect spillage. Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances.

STORAGE INSTRUCTIONS / SHELF LIFE

Protect from exposure to light. May polymerize spontaneously. Store between 4°C and 25°C (39°F-77°F). Close after use. Do not use after the expiration date. See imprint/label for expiration date.

DISPOSAL INSTRUCTIONS

Dispose of contents/container in accordance with local regulations.

IMPORTANT PRODUCT INFORMATION

R_x Federal (U.S.) law restricts the finished device to sale by or on the order of a dentist.

The product was developed for use in dentistry and must be used according to the instructions for use. Serious incidents must be reported to DeltaMed GmbH and the responsible authority.

- LOT** Batch code
- REF** Catalogue number
- MD** Medical device
- Use-by date
- Consult instructions for use
- Temperature limit
- Keep away from sunlight
- Stir before use
- Manufacturer
- Manufacturing Date
- Distributor
- Exclamation mark
- Environmental hazard

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Rapid Shape GmbH