Lucitone Digital Print Denture™ System

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This document contains additional instructions for use information. Refer to the Lucitone Digital Print Denture System Instructions for Use for complete system information, www.dentsplysirona.com.

en

Indications for Use

- Lucitone Digital Fuse™ Step 1 - 3D Tooth Conditioning Agent is indicated for use in enhancing the bond of denture teeth to denture base and denture base to denture base.

Warnings / Adverse Reactions

- Lucitone Digital Print 3D Denture Base and Lucitone Digital Value 3D Economy Tooth & Trial Placement are suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin is suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 3 - 3D Sealer (liquid/vapor) is Toxic to aquatic life with long-lasting effects.
- Lucitone Digital Fuse Step 3 - Total 3D Sealer (liquid/vapor) May damage fertility or the unborn child.

STEP-BY-STEP INSTRUCTIONS

CAUTION:

- Validated printers for use with the Lucitone Digital Print Denture System: Carbon® M-Series, Asiga MAX™ UV, and Asiga PRO 4K™. [1,2]
- Validated post-processing units: Dentsply Sirona inLab Speedcure® Processing Unit (REF #9494800 and 9494801), Dentsply Sirona Digital Cure - Large Capacity Unit (REF #909170 and 909180) and Dentsply Sirona Digital Cure (REF #909188).

A. PRINT

Note: Ensure the latest software version load for complete material availability.

B. RECYCLE

Table I

Printer	Material	Recycle Timeframe
Asiga MAX UV & PRO 4K	Lucitone Digital IPN - 3D Premium Tooth	up to 3 months

E. CURE

CAUTION: Only use validated post-processing units: Dentsply Sirona inLab Speedcure® Processing Unit (Option 1), Dentsply Sirona Digital Cure - Large Capacity Unit (Option 2) or Dentsply Sirona Digital Cure (Option 3). Read the curing unit User Manual prior to operating the unit.

Cure UNIT Option 3: Dentsply Sirona Digital Cure

CAUTION:

- Do not stack appliances - correct positioning is a requirement for proper cure.
- Appliance quantity per cure cycle varies by appliance sizes - do not exceed five (5) appliances per cure cycle. Regardless of the appliance quantity, always place appliances with the anterior facing the center of the turntable.
- Allow appliances cleaned with Isopropyl alcohol to dry completely before post-processing due to flammability.
- Once the program cycle starts do not stop and restart, open the door, or add appliances; a complete, uninterrupted program run is required.
- Regardless of the appliance quantity, always use the Lucitone menu selection (preset). Do not use the custom menu selection to process Lucitone Digital Print Denture appliances.
- Failure to follow specific cure instructions will result in poor appliance strength, bond, and accuracy.

Table IX

Menu Selection	Total Time
Lucitone	Approximately 60 minutes (maximum 70 minutes)
1	Place the appliance(s) on the turntable, occlusal side up, and the anterior facing the center of the turntable. Load multiple appliances using the Start/stop button to the right of the menu dial. Hold the power button to rotate the turntable and place the appliances. Ensure the appliance(s) do not overlap or extend beyond the turntable. Do not load more than five appliances per cycle. IMPORTANT: When curing one appliance, do not place the appliance in the center of the turntable covering the turntable spindle. Ensure the appliance is facing anterior towards the center and approximately 1 cm from the turntable spindle/center.
2	In the home menu screen, use the menu dial to toggle and select (push menu dial) the "Preset Cycle" setting and then the "Lucitone" program. Hold the Start/stop button until the unit beeps to start the program. The program will automatically run and advance through the following steps: LED Test Cycle - 20 seconds, Tack curing - 12 seconds, Heating - 10 to 15 minutes, Curing - 25 minutes, and Cooling - up to 30 minutes. Synchronous dots will appear behind the step name on the menu screen indicating the program operation. Note: the chamber is lit during the Tack Curing and Curing steps only. Following the final step, the menu screen "Program Complete" will appear indicating the program completion. The appliances are ready for removal.

G. REPAIRS

Table IX

	Surface	Fracture*	Tooth - IPN 3D
8	IPN 3D: Apply a thin layer of Lucitone Digital Fuse Step 3 (3D Sealer or Total) to the margins and interproximals using a flocked applicator tip. Surface & Fracture: Apply Lucitone Digital Fuse Step 3 (3D Sealer or Total) to the repair area.		
9	inLab Speedcure: Place the appliance repair side up in the unit. Select the menu option "Lucitone Digital" and then "Repair." Program cycle: 10 minutes, then 3-minute cooldown. Caution: Position the repair surface up on the turntable and only cure one denture at a time. DS Digital Cure - Large Capacity: Place the appliance repair side up in the unit. Select the "Lucitone Digital Appliances" menu selection. The program cycle will automatically run for 90 minutes - "HEATING" will display on the screen when the cycle starts. Note: The unit may require warming prior to the cycle start - "PRE-HEATING" will display on the screen until the unit reaches temperature. Cooldown Cycle: Wait until the temperature indicator in the upper right corner of the home screen reaches 50°C before removing and finishing the appliance(s).		
10	Finish and polish the denture using conventional techniques.		

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