Intended Use

Indication: RS Model Tan is a material based on acrylic esters for manufacturing of all types of dental models within the 3D printers of series D II and D+ of Rapid Shape.

Indication for Use: RS Model Tan is a light curing resin intended for manufacturing of dental models.

Contraindication: RS Model Tan is contraindicated for individuals with known allergy to acrylic esters and acrylates.

Patient Target Group: Patients for whom dental treatment is intended.

Intended Users: The dental technicians or specialist staff in the dental laboratory is responsible for the production of dental models in 3D printing.

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Software: Autodesk Netfabb®

Post-Processing Units: RS wash + RS cure

Basic Material Properties

Characteristics/Properties	Units	Specifications
Before curing (liquid state)		
Viscosity	cps	450-500 (at 25°C)
Density	g/cm ³	1.135
After curing		
Tensile strength	MPa	40-45
Elongation at break	%	8-15
Izod impact	J/m	55-70
Hardness	D	80-85
Glass transition	°C	70-85

Specific Manufacturing Considerations

1. Digital Dental Model File

1.1 File Format: STL file

1.2 Digital Design: dental model

2. Rapid Shape Printers and Software

2.1 Hardware

Please refer to the printer's manual for this information.

2.2 Rapid Shape Printer Software (Autodesk Netfabb®)

Please refer to the printer's manual for this information.

2.3 Environmental Conditions

a) Print Temperature: Please refer to the printer's manual for this information.

b) Humidity: ≤ 50 %

3. Notification

The specifications have been validated using printer software Autodesk Netfabb® and aforementioned printers.

Processing

General Notes

- Make sure that you work as clean as possible. Dirty material reservoirs or machines can cause deformations and/or failed print jobs.

- To avoid damage of material quality, do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or color variations of the RS Model Tan.

Preparation

- During the application of RS Model Tan, we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc.
- Before using the RS Model Tan material, mix it well in its original packaging and each subsequent time you use it. Shake well for a few minutes or use an ultrasonic equipment to thoroughly mix the material. Color deviations and/or print failures may occur when mixed insufficiently.
- Carefully pour the material into the pre-determined reservoir of the printer.
- Remove all bubbles with a cleaned object.
- Please see further steps within the printer's manual.
- Processing temperature: 25±3°C
- Select the buildstyle for RS Model Tan in the printer software. Make sure that all software is up to date.

Further Processing

- After completion of the building process, direct further processing (see following steps) is recommended.
- After the platform has been raised, a dripping time of approx. 10 minutes is recommended.
- Remove the build platform from the printer and place it on a paper towel with the build parts facing upwards. Remove the printed parts with a suitable knife or spatula.

Cleaning/Washing

- Place the build part in the RS wash unit and select the respective parameter set for RS Model Tan. Start the cleaning cycle. (Alternatively, rinse the build part for 3 minutes in 97 % isopropanol (IPA) in an ultrasonic bath to remove any excess material. Dry the cleaned build part with compressed air afterwards.)

Post Curing

- The post curing process is an essential step to ensure a full polymer conversion and to obtain the required mechanical properties of the end-product.
- Option 1: Place the cleaned build part in the RS cure unit and select the respective parameter set for RS Model Tan. Start the curing cycle.
- Option 2: Place the cleaned build part in the Otofash G171 (NK Optik) and cure the part with 2000 flashes each for top and bottom of the object (Σ4000) directly on the grating without insert shell.

Others

- Remove the reservoir from the printer and filter the resin through a fine 190 Micron paint strainer, if the print has failed partially or completely, or particulates of polymerized residues are visible in the reservoir or stick to the bottom.
- To avoid bubbles, let the material rest for 1 to 2 hours prior to use, after filtration.
- Discard and replace the RS Model Tan material with a new batch, if contamination, evident gelation or polymerization is observed after filtering.
- Do not mix different lot numbers of RS Model Tan.

Caution

The specifications have been validated using printer software Autodesk Netfabb® and aforementioned printers and indicated process parameters, herein. Any unauthorized changes to the process equipment, parameters or software may result in a device that is out of specifications.

Disposal

Dispose RS Model Tan in accordance with local/regional/national/international regulations. See also information in safety data sheet.

Hazards Identification

Before use observe the Material Safety Data Sheet of RS Model Tan.

Hazard Statements: H315 Causes skin irritation. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P280 Wear protective gloves/eye protection/face protection. P273 Avoid release to the environment. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P333+P313 If skin irritation or rash occurs: Get medical advice/ attention. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: Acrylic esters (acrylates); photoinitiators; pigments; additives.

Lot Number / Expiry Date

The lot number and the expiry date are indicated on each RS Model Tan packaging. In case of complaints, please always indicate the lot number of the product. Do not use the product after the expiry date.

Storage Conditions

Store the RS Model Tan in its original packaging at room temperature (do not exceed 25°C) and in a dry and dark area, protected from light and moisture.

Forms of Delivery

REF RS001781, RS Model Tan, 1000 g

REF RS001824 RS Model Tan, 10 kg

Glossary

	Manufacturer
	Batch Code
	Item Number
	Protect from Sunlight
	Temperature Limitation
	Expiry Date
	Please note Instructions for Use

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Rapid Shape GmbH is not liable for any damages caused by improper application of the material.

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Rapid Shape GmbH n'est pas responsable des dommages causés par une application incorrecte du matériau.