Lucitone Digital Print Denture™ System for Primeprint

3D Printing Materials
3D-Druckmaterialien
Matériaux pour impression en 3D
Materiali per la stampa 3D
Materiales para la impresión en 3D
Materiais de Impressão 3D
3D-printmaterialen
3D-utskriftsmaterial
3D-printmaterialer
Materialer for 3D-utskrift
Materiály pro 3D tisk

Instructions for Use.....	English
Mode d'emploi.....	French
Istruzioni per l'uso	Italian
Gebrauchsanweisung	German
Instrucciones de uso.....	Spanish
Instruções de utilização.....	Portugese
Gebruiksaanwijzing	Dutch
Bruksanvisning.....	Swedish
Brugsanvisning	Danish
Bruksanvisning.....	Norwegian
Návod k použití.....	Czech

Indications for Use

- Lucitone Digital Print™ 3D Denture Base is a light-cured resin indicated for the fabrication of denture bases in dental laboratories, including full and partial dentures* and implant overdentures.
- Lucitone Digital IPN™ 3D Premium Tooth is a light-cured resin intended for printing denture teeth.
- Lucitone Digital Value™ 3D Economy Tooth & Trial Placement is used as a try-in material for evaluation prior to fabrication of the final restoration and is intended for printing full arch and tooth segments.
- Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent is indicated for use in enhancing the bond of denture teeth to denture base and denture base to denture base.
- Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin is utilized as an aid in bonding denture teeth to denture base as well as repair using traditional techniques.
- Lucitone Digital Fuse™ Step 3 – 3D Sealer is a light-cured sealant that produces a smooth, glossy surface finish on the denture.
- Lucitone Digital Fuse™ Step 3 – Total 3D Sealer is a light-cured sealant that produces a smooth, glossy surface finish on the denture.

*Partial and full dentures are replacement for patients with missing teeth.†

†Statement added for EU MDR alignment.

Performance

- Meets requirements of ISO 20795-1, Type 4

Intended Purpose and Benefit

- Device is intended for fabrication of removable dental prostheses to restore function and aesthetics.
- Patient Target group(s) - Patients with missing teeth

Contraindications

- Lucitone Digital Print Denture System components are contraindicated for patients with known hypersensitivity or severe allergic reactions to any of the components.

Warnings

- Lucitone Digital Print Denture System materials contain components that may cause skin dryness or irritation, sensitization (allergic contact dermatitis), or other allergic reactions in susceptible persons.
 - Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
 - Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
- Lucitone Digital Print Denture System materials: Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea, and irritation of the eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects.
 - Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
 - Ingestion: Contact your regional Poison Control Center immediately.

- Lucitone Digital Print 3D Denture Base and Lucitone Digital Value 3D Economy Tooth & Trial Placement are suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent liquid is highly flammable. Keep container tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization, or other allergic reactions.
- Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin is suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 3 - 3D Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. Suspected of damaging fertility or the unborn child. Toxic to aquatic life with long-lasting effects. Avoid release into the environment.
- Lucitone Digital Fuse Step 3 - Total 3D Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. May damage fertility or the unborn child. Harmful to aquatic life with long-lasting effects. Avoid release into the environment.
- Dispose of contents and containers in accordance with local and national regulations.

Precautions

- These products are intended to be used only as specifically outlined in the Instructions for Use (IFU). Any use of these products inconsistent with the IFU is at the discretion and sole responsibility of the practitioner.
- Store resins and system components at room temperature 60°F - 80°F (16°C - 27°C) and avoid direct sunlight. Keep containers closed when not in use. Do not use products after expiration dates. Dispose of contents and containers in accordance with local and national regulations.
- Lucitone Digital Print 3D Denture Base, Lucitone Digital IPN 3D Premium Tooth, Lucitone Digital Print 3D Economy Trial Placement, and Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin: When washing a 3D printed appliance with a solvent or polishing the appliance, use in a properly ventilated environment with proper protective masks and gloves. Dispose of unused resins in accordance with local and national regulations.
- Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent: Use in a ventilated workspace. Handle teeth with solvent-resistant gloves and tweezers. The product is sensitive to light - minimize exposure to light and keep the lid on the container to maintain its effective use. The product is very volatile and flammable. Keep away from heat, sparks, and flame. Do not add Lucitone Digital Fuse Step 1 to a warm metal container or while the container is on the electric hotplate.
- Lucitone Digital Fuse Step 3 - 3D Sealer and Lucitone Digital Fuse Step 3 - Total 3D Sealer: Use in a ventilated workspace. Keep bottle tightly closed when not in use. Lucitone Digital Fuse Step 3 sealers are volatile and highly flammable. Keep away from heat, sparks, and flames. The product is sensitive to light - minimize exposure to light and keep the lid on the container to maintain its effective use.

Adverse Reactions

- Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
- Proper ventilation and personal protective equipment should be used when grinding devices as the particulate generated during grinding may cause respiratory, skin and/or eye irritation.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local and national regulations.

Requirements for Denture Fabrication

- File(s) in .stl format
- Lucitone Digital Print 3D Denture Base - cartridge for Primeprint Solution Available in shades: Original, Light, Light Reddish Pink, and Dark Reddish Pink
- Lucitone Digital IPN 3D Premium Tooth - cartridge for Primeprint Solution Available in shades: A1, A2, A3, A3.5, B1, B2, C1, and BL3
- Lucitone Digital Print 3D Economy Trial Placement - cartridge for Primeprint Solution Available in shade: A1
- Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin
- Lucitone Digital Fuse Step 3 - Total 3D Sealer
- Primeprint Box
- Material Unit(s)
- $\geq 99\%$ Isopropyl Alcohol (IPA)
- DEG monobutyl ether solvent (ex. InovaPrint Wash)
- UV Tack-Cure Light
- Biobond Brush, #4 Flat
- 4x4 Gauze Pads
- Denture holder
- Individual Curing Carrier
- InLab CAM software - version 2.24 or higher with CAM menu availability for the Lucitone Digital Denture System
- Primeprint
- Primeprint PPU (Post Processing Unit)

Note: Primeprint cartridges are not refillable

STEP-BY-STEP INSTRUCTIONS

CAUTION:

- Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification and is not recommended. Contact Dentsply Sirona for a list of compatible components.
 - Lucitone Digital Print Denture System workflow validations include processes, materials, equipment, CAD/CAM parameters, and software. Deviations may result in a device that is out of specification. Contact Dentsply Sirona for a list of compatible components.
 - Instructions in this IFU feature processes related to the Primeprint printing and processing solution, and feature printed dentures and printed denture teeth. For alternative printer systems and denture workflows, repairs and relines visit dentsplysirona.com - [Lucitone Digital Print Denture Curriculum](#)
 - For detailed step-by-step instructions with images, refer to the Dentsply Sirona Print to Finish Illustrated Technique Guide for Primeprint.
 - Validated printers and processing units covered by these instructions for the Lucitone Digital Print Denture System:
 - Primeprint
 - Primeprint PPU (Post Processing Unit)
 - Use proper ventilation. Wear protective gloves, eye protection, and protective clothing.
-

A. PRINT

CAUTION:

- Instructions are applicable to Lucitone Digital Print Denture System materials only.
 - Design Requirement Maxillary denture/try-in requires 3 T-bars and an anterior extension bar directly behind centrals. Mandibular denture/try-in requires 1 T-bar and 2 support bars. See the Digital Denture Design in 3Shape Illustrated Technique Guide for instructions.
 - Refer to the Primeprint Solution Instructions for Use for complete printer information.
-

1. Use the inLab CAM software to upload the .stl file(s) for the print job. Select “+” to create a “New Job”.
2. Select Primeprint and the Machine Type, select the Item Type to be printed, Manufacturer, Material, and Production Method
***Note:** Ensure the latest software version is loaded for complete material availability.
Options for Lucitone Digital Denture System include:
 - Denture Base: Primeprint Lucitone Denture Base
 - Try-In: Primeprint Lucitone Value Try-In
 - Teeth: Primeprint Lucitone IPN Teeth
3. Mark the areas of the design where no support structures are allowed in the “Add Additional Information” tab.
4. Select the fast forward button for preparation of the print file. If auto-orient feature does not position bases with posterior facing build platform, then manual orientation may be required. Please see the Illustrated technique guide for Manual preparation.

5. Load Primeprint with a material unit that contains the selected resin cartridge and a Primeprint Box.
 6. Select: Start Production
- * For detailed print instructions refer to the Dentsply Sirona Illustrated Technique Guide for Primeprint

B. CLEAN

CAUTION:

- Use proper ventilation. Wear protective gloves, eye protection, and protective clothing.
 - When manually cleaning the residual resin with an 99% IPA-soaked brush, do not expose the part for extended periods of time in IPA. Quickly clean the part surfaces and let air dry. Prolonged Isopropyl Alcohol (IPA) exposure may change appliance physical properties.
-

1. Remove Primeprint Box from Primeprint and put it into Primeprint PPU.
2. Load Primeprint PPU with the washing containers required for the printed resin and unlock their lids as soon as they are inside the device. One washing container is filled with DEG monobutyl ether solvent (ex. InovaPrint Wash) and one is filled with 99% Isopropanol (IPA).
3. Start washing
4. When washing cycle is completed, take the Primeprint Box out of the Primeprint PPU.
5. Remove the print job from the building platform. Detach the parts by hand, using moderate pressure, or use a spatula/scrapper.
6. Remove the support structures using hand pressure or a cutting tool. (Do not remove the bar supports (Denture Bases/Try-ins) - removal will occur after the post-cure step)
7. Use a brush soaked with clean 99% IPA to remove any residual resin.
8. Use compressed air to thoroughly dry part(s) before the Fuse process.

C. FUSE (BOND)

For Trial Placement (Try-In Appliances) skip to section D.

CAUTION:

- Use proper ventilation. Wear protective gloves, eye protection, and protective clothing.
-
1. Prepare printed denture teeth and printed denture by removing the remaining support structures. A carbide bur or cutting disc is recommended for removal of the remaining high spots of the support structures.
 2. Use a diamond disc to open the interproximal areas of the denture teeth. Shape and contour the teeth and denture as desired.
 3. Remove any remaining resin using a small brush and 99% IPA.
 4. Confirm the fit of the denture teeth in the denture base.
Note: Some gaps may be visible within the pockets - these gaps will be filled later with Lucitone Digital Fuse Step 2 resin.
 5. Apply Lucitone Digital Fuse Step 2 into the denture pockets:

For tooth segment:

- 1) Apply Lucitone Digital Fuse Step 2 resin to the tooth quadrant pocket on the denture base.
- 2) Insert a tooth segment, fully seat, and clean excess resin with a gauze.
- 3) Tack cure the buccal and lingual surfaces with the UV Tack-Cure Light for ten seconds, two teeth at a time.
- 4) Repeat the above steps for the remaining tooth segments.

For full tooth arches:

- 1) Apply Lucitone Digital Fuse Step 2 resin to the full arch indentation on the denture base.
 - 2) Insert the full arch, fully seat, and clean excess resin with gauze.
 - 3) Tack cure the buccal and lingual surfaces with the UV Tack-Cure Light for ten seconds, two teeth at a time.
6. Apply a thin layer of Lucitone Digital Fuse Step 3 Total sealer to the entire tooth structures, including the cervical areas, using a flat brush.
 7. The appliance is now ready to cure.

Optional (not required): Apply the Lucitone Digital Fuse Step 3 Total sealer to all or part of the denture base surface (including the intaglio/tissue) per user preference.

D. CURE - Using Primeprint PPU

CAUTION:

- Read the Primeprint PPU curing unit User Manual prior to operating the unit.
 - Failure to follow specific cure instructions will result in poor appliance strength, bond, and accuracy.
-
1. After denture base and denture teeth have been fused, place one denture on one denture holder, with another support bar leaning against it (see the ITG for visual guidance.)
 2. Place a maximum of two arches with denture holders in the Individual Curing Carrier.
 3. Load Primeprint PPU with the Individual Curing Carrier and open the nitrogen bottle.
 4. Close the door of the Primeprint PPU.
 5. Select “Lucitone Digital Fuse” to run the tack-cure cycle.
 6. Select appropriate matching job to begin final cure cycle.
 7. Allow appliance(s) to fully cool prior to handling.

E. FINISH

CAUTION:

- Use proper ventilation and personal protective equipment when grinding devices as the particulate generated during grinding may cause respiratory, skin, and/or eye irritation.
 - Do not steam clean Lucitone Digital Print Denture appliances.
-

1. Remove T-bar and bar supports using a cut-off disc or similar tool. Grind away any remaining support structure.
2. Finish, festoon, stipple, and polish per laboratory preference.

IMPORTANT: Avoid extensive polishing for dentures designed with the Dentsply Sirona Highly Characterized Denture Tooth Libraries to preserve the premium features and characterization.

F. RELINE Hard Reline

Refer to the Lucitone® HIPA - High Impact Pour Acrylic Instructions for Use "Indirect Reline Procedure" section for step-by-step instructions.

Soft Reline

The following materials are compatible for use with Lucitone Digital Print Dentures:

- Tokuyama SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium [1]
- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft, and GC RELINE™ II Extra Extra Soft [2]

QUESTIONS AND SUPPORT

USA - Email: [dsgdigitaldentures@dentsplysirona.com](mailto:digitaldentures@dentsplysirona.com) or Phone: 800-243-1942.

CANADA - Contact your local Dentsply Sirona Representative.

OUTSIDE OF NORTH AMERICA - Contact your local Dentsply Sirona Representative.

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- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft a GC RELINE™ II Extra Extra Soft [2]