Lucitone Digital Print Denture™ System

3D Printing Materials

3D-Druckmaterialien

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Materiali per la stampa 3D

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3D-printmaterialen

3D-utskriftsmaterial

3D-printmaterialer

Materialer for 3D-utskrift

Materiály pro 3D tisk

Materiale de imprimare 3D

Instructions for Use.....	English
Mode d'emploi	French
Istruzioni per l'uso.....	Italian
Gebrauchsanweisung.....	German
Instrucciones de uso.....	Spanish
Instruções de utilização	Portuguese
Gebruiksaanwijzing.....	Dutch
Bruksanvisning.....	Swedish
Brugsanvisning	Danish
Bruksanvisning.....	Norwegian
Návod k použití.....	Czech
Instrucțiuni de utilizare.....	Romanian
Návod na použitie	Slovakian

Lucitone Digital Print Denture™ System

3D Printing Materials

3D-Druckmaterialien

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Materiali per la stampa 3D

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Instructions for Use

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Mode d'emploi

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Indications for Use

- Lucitone Digital Print™ 3D Denture Base is a light-cured resin indicated for the fabrication of denture bases including full and partial dentures* and implant overdentures.
- Lucitone Digital IPN™ 3D Premium Tooth is a light-cured resin intended for printing denture teeth.
- Lucitone Digital Value™ 3D Economy Tooth & Trial Placement is used as a try-in material for evaluation prior to fabrication of the final restoration or a temporary denture.
- Lucitone Digital Value™ 3D Economy Tooth & Trial Placement is intended for printing full arch and tooth segments.
- Lucitone Digital Fuse™ Step 1 – 3D Tooth Conditioning Agent is indicated for use in enhancing the bond of denture teeth to denture base and denture base to denture base.
- Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin is utilized as an aid in bonding denture teeth to denture base as well as repair using traditional techniques.
- Lucitone Digital Fuse™ Step 3 – Total 3D Sealer is a light-cured sealant that produces a smooth, glossy surface finish on the denture.

*Partial and full dentures are replacement for patients with missing teeth.†

†Statement added for EU MDR alignment.

Performance

- Meets requirements of ISO 20795-1, Type 4

Intended Purpose and Benefit

- Device is intended for fabrication of removable dental prostheses to restore function and aesthetics.
- Patient Target group(s) - Patients with missing teeth

Contraindications

- Lucitone Digital Print Denture System components are contraindicated for patients with known hypersensitivity or severe allergic reactions to any of the components.

Warnings

- Lucitone Digital Print Denture System materials contain components that may cause skin dryness or irritation, sensitization (allergic contact dermatitis), or other allergic reactions in susceptible persons.
 - Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
 - Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.
- Lucitone Digital Print Denture System materials: Avoid inhalation or ingestion. High vapor concentration might induce headache, nausea, and irritation of the eyes and respiratory system. Work in a well-ventilated area. Avoid contact with eyes. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects.
 - Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.
 - Ingestion: Contact your regional Poison Control Center immediately.

- Lucitone Digital Print 3D Denture Base and Lucitone Digital Value 3D Economy Tooth & Trial Placement are suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent liquid is highly flammable. Keep container tightly closed and use/store away from sources of ignition. Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization, or other allergic reactions.
- Lucitone Digital Fuse Step 2 – 3D Denture Bonding Resin is suspected of damaging fertility or the unborn child.
- Lucitone Digital Fuse Step 3 - Total 3D Sealer (liquid/vapor) is highly flammable. Keep containers tightly closed and use/store away from heat or sources of ignition. May damage fertility or the unborn child. Harmful to aquatic life with long-lasting effects. Avoid release into the environment.
- Dispose of contents and containers in accordance with local and national regulations.

Precautions

- These products are intended to be used only as specifically outlined in the Instructions for Use (IFU). Any use of these products inconsistent with the IFU is at the discretion and sole responsibility of the practitioner.
- Store resins and system components at room temperature 60°F - 80°F (16°C – 27°C) and avoid direct sunlight. Keep containers closed when not in use. Do not use products after expiration dates. Dispose of contents and containers in accordance with local and national regulations.
- Lucitone Digital Print 3D Denture Base, Lucitone Digital IPN 3D Premium Tooth, Lucitone Digital Print 3D Economy Tooth & Trial Placement, and Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin: When washing a 3D printed appliance with a solvent or polishing the appliance, use in a properly ventilated environment with proper protective masks and gloves. Dispose of unused resins in accordance with local and national regulations.
- Lucitone Digital Fuse Step 1 – 3D Tooth Conditioning Agent: Use in a ventilated workspace. Handle teeth with solvent-resistant gloves and tweezers. The product is sensitive to light - minimize exposure to light and keep the lid on the container to maintain its effective use. The product is very volatile and flammable. Keep away from heat, sparks, and flame. Do not add Lucitone Digital Fuse Step 1 to a warm metal container or while the container is on the electric hotplate.
- Lucitone Digital Fuse Step 3 - Total 3D Sealer: Use in a ventilated workspace. Keep bottle tightly closed when not in use. Lucitone Digital Fuse Step 3 Total 3D Sealer is volatile and highly flammable. Keep away from heat, sparks, and flames. The product is sensitive to light - minimize exposure to light and keep the lid on the container to maintain its effective use.

Adverse Reactions

- Allergic contact dermatitis and other allergic reaction may occur in susceptible individuals. Product may cause skin dryness, sensitization or other allergic reactions.
- Proper ventilation and personal protective equipment should be used when grinding devices as the particulate generated during grinding may cause

respiratory, skin and/or eye irritation.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local and national regulations.

STEP-BY-STEP INSTRUCTIONS

CAUTION:

- Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification and is not recommended. Contact Dentsply Sirona for a list of compatible components.
 - Lucitone Digital Print Denture System workflow validations include processes, materials, equipment, CAD/CAM parameters, and software. Deviations may result in a device that is out of specification. Contact Dentsply Sirona for a list of compatible components.
 - Validated printers for use with the Lucitone Digital Print Denture System: Carbon® M-Series, Asiga MAX™ UV, Asiga PRO 4K™, Asiga Ultra, HeyGears® UltraCraft A2D and HeyGears® UltraCraft A2D 4K [1,2,7].
 - Validated hotplates: VWR® Hotplates (REF #906235 and 906237) and OHAUS® Guardian™ 5000 Hotplates (REF #906236 and 906238). [3,4]
 - Validated post-processing units: Dentsply Sirona inLab Speedcure® Processing Unit (REF #9494800 and 9494801), Dentsply Sirona Digital Cure - Large Capacity Unit (REF #909170 and 909180) and Dentsply Sirona Digital Cure (REF #909188).
 - Use proper ventilation. Wear protective gloves, eye protection, and protective clothing.
 - For detailed step-by-step instructions with images, refer to the Dentsply Sirona Print to Finish Illustrated Technique Guide for the specific printer manufacturer.
-

A. PRINT

1. Follow the CAD manufacturer's instructions to create a .stl file per design. Use the printer's CAM software to upload the .stl file(s) for the print job.*
Note: Ensure the latest software version load for complete material availability.
 2. Add resin to the printer and start the printer.*
- *For detailed print instructions refer to the printer manufacturer's Instructions for Use/ User Guide and the Dentsply Sirona Illustrated Technique Guide for the specific printer manufacturer.

B. RECYCLE

CAUTION:

- Follow the printer manufacturer's Instructions for Use.
- Do not leave resin in the cassette/vat for an extended period of time; minimize the exposure to ambient light with resin in the cassette/vat, and pour resin back into the bottle immediately after the print job is complete.
- Do not mix resin shades.

1. Filter left-over resin into the original container. Refer to Table I for recycle timeframes.

Table I

Printer	Material	Recycle Timeframe
Carbon M-Series	Lucitone Digital Print - 3D Denture Base	up to 5 months
	Lucitone Digital Value - 3D Economy Tooth & Trial Placement	
	Lucitone Digital IPN - 3D Premium Tooth	
Asiga MAX UV/ PRO 4K/ Ultra/ HeyGears® UltraCraft A2D and UltraCraft A2D 4K	Lucitone Digital Print - 3D Denture Base	up to 3 months
	Lucitone Digital Value - 3D Economy Tooth & Trial Placement	
	Lucitone Digital IPN - 3D Premium Tooth	

C. CLEAN

CAUTION:

- Isopropyl Alcohol is highly flammable. Keep containers tightly closed and use/store away from sources of ignition.

- Remove the print job from the printer.
- Detach the part(s) by hand, using moderate pressure, or use a spatula/scrapper.
- Break off the part supports. Leave the bar supports in place - removal will occur after the post-cure step.
- Place the part(s) in a plastic container, add $\geq 99\%$ Isopropyl Alcohol (IPA) covering the part(s) and seal with the container lid. **Note:** Use new IPA only for both cycles.
Cycle 1: Place the container in an ultrasonic bath for two (2) minutes. Remove the container and clean the part(s) using a brush soaked with $\geq 99\%$ IPA to remove the residual resin.
Cycle 2: Using new IPA, return the container to the ultrasonic bath for one (1) minute. Once the cycle is complete, remove the container and check the part(s) as it is critical to ensure there is no residual resin. Use a brush soaked with $\geq 99\%$ IPA to remove any residual resin.
- Use compressed air to thoroughly dry all parts before the Fuse step.

D. FUSE (BOND)

For Trial Placement (Try-In Appliances) skip to section E.

Table II

Summary of Fuse Steps	Digital Tooth Materials			
	Lucitone Digital IPN	Lucitone Digital Value	DS Multilayer PMMA Disc	IPN 3D
Tooth Preparation for Milled Denture Teeth	Not Applicable	Not Applicable	Remove the milling supports, open the interproximals, shape/contour the teeth, and clean in a ultrasonic bath.	Not Applicable
Tooth Preparation for Printed Denture Teeth	Remove the printing supports, open the interproximals, shape/contour the teeth, and clean debris with IPA.	Remove the printing supports, open the interproximals, shape/contour the teeth, and clean debris with IPA.	Not Applicable	Not Applicable

Dry Fit	Place the tooth segments, full arch, or single teeth into the denture base. Some gaps may be visible within the pockets - these gaps will be filled with Fuse Step 2 resin.	Place the tooth segments or full arch into the denture base. Some gaps may be visible within the pockets - these gaps will be filled with Fuse Step 2 resin.	Place the tooth segments or full arch into the denture base. Some gaps may be visible within the pockets - these gaps will be filled with Fuse Step 2 resin.	Place each tooth into the respective denture base pocket to verify fit.
Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent	Not Applicable	Not Applicable	See Table III DS Multilayer PMMA	See Table III IPN 3D
Lucitone Digital Fuse Step 2 - 3D Denture Bonding Resin	See Table IV Single Teeth, Segments or Full Arch	See Table IV Segments or Full Arch	See Table IV Segments or Full Arch	See Table IV Single Teeth
Lucitone Digital Fuse Step 3 Total 3D Sealer	See Table V Lucitone Digital IPN	See Table V Lucitone Digital Value	See Table V DS Multilayer PMMA/IPN 3D	See Table V DS Multilayer PMMA/IPN 3D

Lucitone Digital Fuse Step 1 - 3D Tooth Conditioning Agent
For printed denture teeth skip to Lucitone Digital Fuse Step 2 - Table IV

CAUTION:

- Use proper ventilation. Wear protective gloves, eye protection, and protective clothing.
- Keep the metal container and the foam insert free of any solid buildup. Buildup can compromise the tooth bond strength.

Table III

	DS Multilayer PMMA	IPN 3D
1	Plug in the hotplate and set the temperature per below. Preheat for 10 minutes. VWR Hotplates (REF # 906235 and 906237): Temperature setting 90 °C OHAUS Hotplates (REF # 906236 and 906238): Temperature setting 60 °C	
2	Shake the bottle and pour Step 1 into the metal container.	

3	Place the tooth segments/arch in the metal container - ridge lap down & complete submersion in liquid. Place the lid on the container.	Place the foam insert in the metal container, use tweezers to place the teeth into the foam holes - ridge lap down & complete submersion in liquid. Place the lid on the container.
4	Place the metal container (inside the metal ring for VWR models only) on the hotplate until the temperature strip reaches 40 °C (Approx. 4-6 minutes).*	
5	Remove the teeth/segments/arch from the metal container using tweezers. Air-dry for 2 minutes. Do not wipe.	
*Failure to reach 40 °C or using higher warming temperatures to accelerate the process will compromise tooth bond strength.		

Lucitone Digital Fuse Step 2 - 3D Tooth Bonding Resin

Table IV

	Single Teeth	Segments	Full Arch
1	Apply Step 2 resin to tooth pocket(s) - no more than two teeth at a time - avoiding adjacent pockets.	Apply Step 2 resin to tooth quadrant indentation on the denture base.	Apply Step 2 resin to full arch indentation on the denture base.
2	Insert tooth/teeth , fully seat, and clean excess resin with 4X4 gauze.	Insert segment , fully seat, and clean excess resin with 4X4 gauze.	Insert full arch , fully seat, and clean excess resin with 4X4 gauze.
3	Tack cure the buccal and lingual surfaces with the UV Tack-Cure Light for ten seconds, two teeth at a time.		
4	Repeat above steps 1-3 for any remaining teeth.	Repeat above steps 1-3 for remaining segments.	Not Applicable
NOTE: Fuse Step 2 is available in all Lucitone Digital Print 3D Denture Base shades. Original shade may be used as a universal shade option except for with dark reddish pink.			

Lucitone Digital Fuse Step 3 Total 3D Sealer

CAUTION: When using the Dentsply Sirona Digital Cure - Large Capacity Unit for post-processing, complete the following step after the sealer application: Tack cure the sealer for one (1) minute using the UV Tack-Cure Light. Rotate the arch ensuring light exposure to all tooth surfaces for a uniform tack cure. The sealer will be slightly sticky/tacky to the touch - avoid excess handling. Following the Cure (post-processing) step, the appliances will be fully cured and smooth to the touch.

Table V

Sealer	Lucitone Digital IPN	Lucitone Digital Value	DS Multilayer PMMA/IPN 3D
Lucitone Digital Fuse Step 3 - Total	Apply a thin layer of Total sealer to the entire tooth structures, including the cervicals, using a flat brush. Optional (not required): Apply the sealer to all or part of the denture base surface (including the intaglio/tissue) per laboratory preference.		Apply a thin layer of Total sealer to the cervicals using a flocked applicator tip. Optional (not required): Apply the sealer to all or part of the denture base surface (including the intaglio/tissue) per laboratory preference.

E. CURE

CAUTION:

- Only use validated post-processing units: Dentsply Sirona inLab Speedcure® Processing Unit (Option 1), Dentsply Sirona Digital Cure - Large Capacity Unit (Option 2) or Dentsply Sirona Digital Cure (Option 3). Read the curing unit User Manual prior to operating the unit.
- Failure to follow specific cure instructions will result in poor appliance strength, bond, and accuracy.

Table VI

The post processing units are validated for the following Lucitone Digital Print Denture System materials and appliances:

Application	Dentsply Sirona Digital Materials	Appliance Type
Denture Base	Lucitone Digital Print 3D Denture Base	Denture
Denture Teeth	Lucitone Digital IPN 3D Premium Tooth	
	DS Multilayer PMMA Disc	
	IPN 3D Digital Denture Teeth	
Denture Teeth/ Try-In	Lucitone Digital Value 3D Economy Tooth & Trial Placement	Denture/Try-In

CURE UNIT OPTION 1: Dentsply Sirona inLab Speedcure Processing Unit

CAUTION:

- Do not stack appliances or position appliances in a way to allow for direct contact with post-processing unit bulbs.
- Do not cure more than two (2) appliances at one time or handle appliances before the 3-minute cooldown.

Table VII

Menu Selections		Total Time
Lucitone Digital > Try-In or Denture Base		26 minutes
Try-In		Denture
1	Place up to two (2) appliance(s), occlusal side up, in the processing unit. For two appliance(s), position anterior teeth towards the center of the turntable.	
2	Select the menu option "Lucitone Digital" and then "Try-In." The program will run for ten (10) minutes, followed by a 3-minute cooldown.	Select the menu option "Lucitone Digital" and then "Denture Base." The program will run for ten (10) minutes, followed by a 3-minute cooldown.
3	Flip the appliance(s) and repeat the same "Try-In" menu option.	Flip the appliance(s) and repeat the same "Denture Base" menu option.

CURE UNIT OPTION 2: Dentsply Sirona Digital Cure - Large Capacity Unit

CAUTION:

- Do not stack appliances - correct positioning is a requirement for proper cure.
- Appliance quantity per cure cycle varies by appliance sizes. Appliances must fit within the defined cure zone and have correct positioning. Regardless of the appliance quantity, only cure the appliance(s) in the defined cure zone.
- Allow appliances cleaned with Isopropyl alcohol to dry completely before post-processing due to flammability.
- Once the program cycle starts: do not stop and restart, open the door, or add appliances; a complete, uninterrupted program run is required.
- Regardless of the appliance quantity, always use the Lucitone Digital Appliances menu selection (preset). Do not use the custom menu selection to process Lucitone Digital Print Denture appliances.

Table VIII

Menu Selections		Total Time
Lucitone Digital Appliances		90 minutes
1	Place the appliance(s) in the cure zone, occlusal side up, and the posterior facing the center of the turntable. Ensure the appliance(s) do not overlap or extend beyond the cure zone. Regardless of the appliance quantity, only cure the appliance(s) in the defined cure zone.	
	Cure zone: The space between 5cm and 12cm from the center of the turntable.	
2	Select the "Lucitone Digital Appliances" menu selection. The program cycle will automatically run for 90 minutes - "HEATING" will display on the screen when the cycle starts. NOTE: The unit may require warming prior to the cycle start - "PRE-HEATING" will display on the screen until the unit reaches temperature.	

3	Cooldown Cycle: Wait until the temperature indicator in the upper right corner of the home screen reaches 50 °C before removing and finishing the appliance(s).
---	---

CURE UNIT OPTION 3: Dentsply Sirona Digital Cure

CAUTION:

- Do not stack appliances - correct positioning is a requirement for proper cure.
- Appliance quantity per cure cycle varies by appliance sizes - do not exceed five (5) appliances per cure cycle. Regardless of the appliance quantity, always place appliances with the anterior facing the center of the turntable.
- Allow appliances cleaned with Isopropyl alcohol to dry completely before post-processing due to flammability.
- Once the program cycle starts do not stop and restart, open the door, or add appliances; a complete, uninterrupted program run is required.
- Regardless of the appliance quantity, always use the Lucitone menu selection (preset). Do not use the custom menu selection to process Lucitone Digital Print Denture appliances.
- Failure to follow specific cure instructions will result in poor appliance strength, bond, and accuracy.

Table IX

Menu Selection		Total Time
Lucitone		Approximately 60 minutes (maximum 70 minutes)
1	<p>Place the appliance(s) on the turntable, occlusal side up, and the anterior facing the center of the turntable. Load multiple appliances using the Start/stop button to the right of the menu dial. Hold the power button to rotate the turntable and place the appliances. Ensure the appliance(s) do not overlap or extend beyond the turntable. Do not load more than five appliances per cycle.</p> <p>IMPORTANT: When curing one appliance, do not place the appliance in the center of the turntable covering the turntable spindle. Ensure the appliance is facing anterior towards the center and approximately 1 cm from the turntable spindle/center.</p>	
2	<p>In the home menu screen, use the menu dial to toggle and select (push menu dial) the “Preset Cycle” setting and then the “Lucitone” program. Hold the Start/stop button until the unit beeps to start the program.</p> <p>The program will automatically run and advance through the following steps: LED Test Cycle - 20 seconds, Tack curing - 12 seconds, Heating - 10 to 15 minutes, Curing - 25 minutes, and Cooling - up to 30 minutes. Synchronous dots will appear behind the step name on the menu screen indicating the program operation. Note: the chamber is lit during the Tack Curing and Curing steps only.</p> <p>Following the final step, the menu screen “Program Complete” will appear indicating the program completion. The appliances are ready for removal.</p>	

F. FINISH

CAUTION:

- Use proper ventilation and personal protective equipment when grinding devices as the particulate generated during grinding may cause respiratory, skin, and/or eye irritation.
- Do not steam clean Lucitone Digital Print Denture appliances.

- Remove bar supports using a cut-off disc or similar tool. Grind away any remaining support structure.
- Finish, festoon, stipple, and polish as desired.

Dentures sealed with Lucitone Digital Fuse Step 3 - *Total* require minimal polishing.

IMPORTANT: Avoid extensive polishing for dentures designed with the Dentsply Sirona Highly Characterized Denture Tooth Libraries to preserve the premium features and characterization.

G. REPAIRS

CAUTION: Do not steam clean Lucitone Digital Print Denture appliances.

Table X

	Surface	Fracture*	Tooth - IPN 3D
1	Not Applicable	Tack fractured denture with adhesive glue or sticky wax. Fabricate a putty model with silicone putty.	Not Applicable
2	Grind and roughen the repair surface with a carbide bur.		
3	Not Applicable	Not Applicable	Use Lucitone Digital Fuse Step 1 to treat the replacement tooth/teeth (See Table III). Next, use sticky wax to position the tooth/teeth.
4	Apply Lucitone Digital Fuse Step 1 to the repair area and air-dry for 1 minute.		
5	Not Applicable	Secure the denture on the putty model with sticky wax or similar type material.	Not Applicable
6	Apply Lucitone Digital Fuse Step 2 to the repair area with a brush. Allow the material to flow into the repair area to avoid air pockets - do not brush.		
7	Tack cure with the UV Tack-Cure Light for 10 seconds. Repeat Lucitone Digital Fuse Step 2 and UV light cure to achieve a desired thickness.		

8	<p>IPN 3D: Apply a thin layer of Lucitone Digital Fuse Step 3 (Total 3D Sealer) to the margins and interproximals using a flocked applicator tip. Surface & Fracture: Apply Lucitone Digital Fuse Step 3 (Total 3D Sealer) to the repair area.</p>
9	<p>inLab Speedcure: Place the appliance repair side up in the unit. Select the menu option “Lucitone Digital” and then “Repair.” Program cycle: 10 minutes, then 3-minute cooldown. Caution: Position the repair surface up on the turntable and only cure one denture at a time.</p> <p>DS Digital Cure - Large Capacity: Place the appliance repair side up in the unit. Select the “Lucitone Digital Appliances” menu selection. The program cycle will automatically run for 90 minutes - “HEATING” will display on the screen when the cycle starts. Note: The unit may require warming prior to the cycle start - “PRE-HEATING” will display on the screen until the unit reaches temperature. Cooldown Cycle: Wait until the temperature indicator in the upper right corner of the home screen reaches 50° C before removing and finishing the appliance(s).</p>
9	<p>DS Digital Cure: In the home menu screen, use the menu dial to toggle and select (push menu dial) the “Preset Cycle” setting and then the “Lucitone” program. Hold the Start/stop button until the unit beeps to start the program.</p> <p>The program will automatically run and advance through the following steps: LED Test Cycle - 20 seconds, Tack curing - 12 seconds, Heating - 10 to 15 minutes, Curing - 25 minutes, and Cooling - up to 30 minutes. Synchronous dots will appear behind the step name on the menu screen indicating the program operation. Note: the chamber is lit during the Tack Curing and Curing steps only.</p> <p>Following the final step, the menu screen “Program Complete” will appear indicating the program completion. The appliances are ready for removal.</p>
10	<p>Finish and polish the denture using conventional techniques.</p>
<p>* Lucitone® HIPA - High Impact Pour Acrylic is an alternative fracture repair technique. Refer to the Lucitone HIPA - High Impact Pour Acrylic Instructions for Use “Repair Procedure” section for step-by-step instructions.</p>	

H. RELINE

Hard Reline

Refer to the Lucitone® HIPA - High Impact Pour Acrylic Instructions for Use “Indirect Reline Procedure” section for step-by-step instructions.

Soft Reline

The following materials are compatible for use with Lucitone Digital Print Dentures:

- Tokuyama SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium [5]
- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft, and GC RELINE™ II Extra Extra Soft [6]

QUESTIONS AND SUPPORT

USA - Email: dsdigitaldentures@dentsplysirona.com or Phone: 800-243-1942.

CANADA - Contact your local Dentsply Sirona Representative.

OUTSIDE OF NORTH AMERICA - Contact your local Dentsply Sirona Representative.

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- Tokuyama SOFRELINER TOUGH® Soft y SOFRELINER TOUGH® Medium [5]
- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft y GC RELINE™ II Extra Extra Soft [6]

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- [1] Carbon is a registered trademark of Carbon, Inc./Carbon ist ein eingetragenes Markenzeichen von Carbon, Inc./Carbon est une marque déposée de Carbon, Inc./Carbon è un marchio registrato di Carbon, Inc./Carbon es una marca registrada de Carbon, Inc.
- [2] Asiga, Asiga MAX, Asiga PRO 4K and Asiga Ultra are trademarks of Asiga./Asiga, Asiga MAX, Asiga PRO 4K et Asiga Ultra sind Markenzeichen von Asiga./Asiga, Asiga MAX, Asiga PRO 4K et Asiga Ultra sont des marques de commerce d'Asiga./Asiga, Asiga MAX, Asiga PRO 4K e Asiga Ultra sono marchi di Asiga./Asiga, Asiga MAX, Asiga PRO 4K y Asiga Ultra son marcas registradas de Asiga.
- [3] VWR is a registered trademark of VWR International, LLC./VWR ist ein eingetragenes Markenzeichen von VWR International, LLC./VWR est une marque déposée de VWR International, LLC./VWR è un marchio registrato di VWR International, LLC./VWR es una marca registrada de VWR International, LLC.
- [4] OHAUS and Guardian are trademarks of OHAUS Corporation./OHAUS und Guardian sind Markenzeichen von OHAUS Corporation./OHAUS et Guardian sont des marques de commerce d'OHAUS Corporation./OHAUS e Guardian sono marchi di OHAUS Corporation./OHAUS y Guardian son marcas registradas de OHAUS Corporation.
- [5] SOFRELINER TOUGH® Soft and SOFRELINER TOUGH® Medium are registered trademarks of Tokuyama Corporation./SOFRELINER TOUGH® Soft und SOFRELINER TOUGH® Medium sind eingetragene Markenzeichen von Tokuyama Corporation./SOFRELINER TOUGH® Soft et SOFRELINER TOUGH® Medium sont des marques déposées de Tokuyama Corporation./SOFRELINER TOUGH® Soft e SOFRELINER TOUGH® Medium sono marchi registrati di Tokuyama Corporation./SOFRELINER TOUGH® Soft y SOFRELINER TOUGH® Medium son marcas registradas de Tokuyama Corporation.
- [6] GC RELINE™ II Soft, RELINE™ II Extra Soft and GC RELINE™ II Extra Extra Soft are trademarks of GC Corporation./GC RELINE™ II Soft, RELINE™ II Extra Soft und GC RELINE™ II Extra Extra Soft sind Markenzeichen GC Corporation./GC RELINE™ II Soft, RELINE™ II Extra Soft et GC RELINE™ II Extra Extra Soft sont des marques de commerce de GC Corporation./GC RELINE™ II Soft, RELINE™ II Extra Soft e GC RELINE™ II Extra Extra Soft sono marchi di GC Corporation./GC RELINE™ II Soft, RELINE™ II Extra Soft y GC RELINE™ II Extra Extra Soft son marcas registradas de GC Corporation.
- [7] HeyGears® UltraCraft A2D and UltraCraft A2D 4K are trademarks of HeyGears®, Inc./HeyGears® UltraCraft A2D und UltraCraft A2D 4K sind Markenzeichen von HeyGears®, Inc./HeyGears® UltraCraft A2D et UltraCraft A2D 4K sont des marques de commerce HeyGears®, Inc./HeyGears® UltraCraft A2D e UltraCraft A2D 4K sono marchi di HeyGears®, Inc./HeyGears® UltraCraft A2D y UltraCraft A2D 4K son marcas registradas de HeyGears®, Inc.

 Manufacturer:
Dentsply LLC
 1301 Smile Way
 York, PA 17404 U.S.A.
 Tel. 800-243-1942
 dentsplysirona.com

EC REP
DeguDent GmbH
 Rodenbacher Chaussee 4
 63457 Hanau-Wolfgang
 Germany
 Tel. +49/6181/5950

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Tabela I

Impressora	Material	Cronograma de reciclagem
Carbon M-Series	Lucitone Digital Print - 3D Denture Base	até 5 meses
	Lucitone Digital Value - 3D Economy Tooth & Trial Placement	
	Lucitone Digital IPN - 3D Premium Tooth	
Asiga MAX UV/PRO 4K/Ultra/HeyGears® UltraCraft A2D e HeyGears® UltraCraft A2D 4K	Lucitone Digital Print - 3D Denture Base	até 3 meses
	Lucitone Digital Value - 3D Economy Tooth & Trial Placement	
	Lucitone Digital IPN - 3D Premium Tooth	

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- GC Corporation GC RELINE™ II Soft, RELINE™ II Extra Soft, a GC RELINE™ II Extra Extra Soft [6]
¹Carbon je registrovaná ochranná známka spoločnosti Carbon, Inc.
³VWR je registrovaná ochranná známka spoločnosti VWR International, LLC.
⁴OHAUS a Guardian sú ochranné známky spoločnosti OHAUS Corporation.
⁵HeyGears® UltraCraft A2D a UltraCraft A2D 4K sú ochranné známky spoločnosti HeyGears®, Inc.