VITA

VITA VIONIC® DENT RESIN

DE	Gebrauchsanweisung Kurzversion	4
EN	Instructions for use Short version see page 120 for US version	8
FR	Mode d'emploi Version courte	12
ES	Instrucciones de uso Versión abreviada	16
IT	Istruzioni d'uso Versione breve	20
BG	Инструкция за употреба Съкратена версия	24
CS	Návod k použití Stručná verze	28
DA	Brugsanvisning kort version	32
EL	Οδηγίες χρήσης Συνοπτικότερη έκδοση	36
ET	Kasutusjuhend lühiversioon	40
FI	Käyttöohje Lyhyt versio	44
HR	Upute za uporabu kratka verzija	48
HU	Használati utasítás Rövid változat	52
JA	使用説明書 ショートバージョン	56
LT	Naudojimo instrukcija Trumpa versija	60

LV	Lietošanas instrukcija Īsā versija	64
NO	Bruksanvisning Kortversjon	68
NL	Gebruiksaanwijzing Korte versie	72
PL	Instrukcja używania Skrócona wersja	76
pt-BR	Instruções de uso Versão resumida	80
pt-PT	Instruções de utilização Versão resumida	84
RO	Instrucțiuni de utilizare Versiunea pe scurt	88
RU	Инструкция по применению Краткая версия	92
SK	Návod na použitie skrátená verzia	96
SL	Navodila za uporabo Kratka različica	100
SV	Bruksanvisning Kort version	104
TR	Kullanım kılavuzu Kısa versiyon	108
UA	Інструкція із застосування Коротка версія	112
ZH	使用說明書 簡短版	116
EN-US	Instructions for use US version Short version	120

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DE

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VITA VIONIC DENT RESIN is a light-cured resin for the fabrication of removable denture teeth using the 3D printing process.

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Dentist, dental technician.

Indication

Permanent denture teeth and long-term provisional crowns and bridges for a wearing time up to 12 months.

Contraindication

Contains (meth)acrylate and phosphine oxides. Some ingredients of VITA VIONIC DENT RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Only insert the product intraorally in completely polymerised state.

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Design process:

1. Prepare the data:
 - CAD design.
Recommended ***bonding gap settings*** in combination with VITA VIONIC BASE RESIN IMPACT: ***150µm***
 - Generation of an STL file
 - Convert the data to the CAM software and create supporting structures. Edit the supporting structures suggested by the software.
 - Conversion to a printer-compatible format.
2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC DENT RESIN from your printer manufacturer.
3. Transfer the prepared data to the 3D printer.

4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
6. Start the printing process.

Post-processing:

7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
9. Remove the outer supporting structure.
10. Remove excess resin with compressed air.

Cleaning:

11. Clean the finished printed objects twice for 1 Min. minutes each time in an ultrasonic bath with fresh isopropanol $\geq 98\%$, according to the device manufacturer's instructions.
12. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning. Store components protected from light at room temperature for 30 minutes before post-exposure.

Bonding:

13. Leave the printed denture base and teeth in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again.
14. Apply liquid VITA VIONIC BASE RESIN IMPACT to the cavities with a small brush or a syringe.
15. Connect the repositioned denture teeth by exposure in a Xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2×1800 flashes under inert gas conditions (nitrogen).

Post-exposure:

16. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The properties and shade effects of the end product depend on a number of factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

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General guidelines

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 50 µm is suitable for the printing process.
- Ideal processing temperature 23°C ± 2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- For long-term provisional crowns and bridges up to 3-units, the minimum occlusal thickness is 1.5 mm, circular 1.5 mm and cervical 1.0 mm. Failure to observe this information may lead to an undesirable result. The temporary crown and bridge restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.

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Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01
After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

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Please read the instructions for use before using the product.

VITA VIONIC DENT RESIN is a light-cured resin for the fabrication of removable denture teeth and temporary restorations such as crowns and bridges. It can be used in combination with a stereolithographic 3D printer using a 385nm light source. The product processing requires a computer-aided manufacturing system (CAD/CAM) that includes the following: scanner, design software, additive printer and post-cure unit. Those devices and software are not part of the product.

CAUTION STATEMENT Caution: Federal Law restricts this device to sale by or on the order of a dentist.

Intended Use

Resin for dental 3D printing.

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users:

Dentist, dental technician (Rx only).

Indication

Permanent denture teeth and long-term provisional crowns and bridges.

Contraindication

Contains (meth)acrylate, phosphine oxides and silanized dental glass. Overall share of inorganic fillers (particle size 1.5 µm) totals 20–40 % by mass. Some ingredients of VITA VIONIC DENT RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Only insert the product intraorally in completely polymerised state.

Product application

Please adhere to the following workflow when using this product:

Design process:

1. Prepare the data:
 - CAD design.
Recommended ***bonding gap settings*** in combination with VITA VIONIC BASE RESIN IMPACT: **150µm**
 - Generation of an STL file
 - Convert the data to the CAM software and create supporting structures. Edit the supporting structures

suggested by the software.

- Conversion to a printer-compatible format.
- 2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC DENT RESIN from your printer manufacturer. Transfer the prepared data to the 3D printer.
- 3. Transfer the prepared data to the 3D printer.
- 4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
- 5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
- 6. Start the printing process.

Post-processing:

- 7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
- 8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
- 9. Remove the outer supporting structure.
- 10. Remove excess resin with compressed air.

Cleaning:

- 11. Clean the finished printed objects twice for 1 Min. minutes each time in an ultrasonic bath with fresh isopropanol $\geq 98\%$, according to the device manufacturer's instructions. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.
- 12. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning. Store components protected from light at room temperature for 30 minutes before post-exposure.

Bonding:

- 13. Leave the printed denture base and teeth in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again.
- 14. Apply liquid VITA VIONIC BASE RESIN IMPACT to the cavities with a small brush or a syringe.
- 15. Connect the repositioned denture teeth by exposure in a Xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2×1800 flashes under inert gas conditions (nitrogen).

Post-exposure:

- 16. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The properties and shade effects of the end product depend on a number of factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

Surface treatment:

17. The denture is now ready for subsequent processing or finishing. Process the surface mechanically (smooth with a milling machine) and polish. Pre-polishing is done using brushes and pre-polishing paste/pumice stone, and high-gloss polishing is done using buffing wheels and high-gloss polishing agents for plastics.

Final cleaning:

18. The cleaning process is carried out by dental laboratories/dental practices in an ultrasonic bath. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.
19. The patient should use lukewarm water, a soft toothbrush, liquid soap or detergent to clean the dentures daily. Rinse well afterwards. Denture cleaning tablets can be used according to the manufacturer's instructions.

Note: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.
Do not use any acids or solvents.

General guidelines

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 50 µm is suitable for the printing process.
- Ideal processing temperature 23°C ± 2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- The 3D printing resin can remain in the vat for a maximum of 14 days (maximum 10 print cycles).
- It is not allowed to pour 3D printing resin from the reservoir of the printer back into the bottle.
- For long-term provisional crowns and bridges up to 3-units, the minimum occlusal thickness is 1.5 mm, circular 1.5 mm and cervical 1.0 mm. Failure to observe this information may lead to an undesirable result. The temporary crown and bridge restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.

Recommended product combinations

For recommended product combinations with the VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Use of this product is only authorized in conjunction with compatible devices. The devices should be set up, validated and maintained in accordance with their labeling and instructions for use. Any unauthorized changes to the process equipment, parameters, or software may result in a product that is out of specification.

An up-to-date list of **validated 3D-printers, washing units and post curing units**, as well as suitable printer templates and post-curing parameters, can be found at www.vita-zahnfabrik.com/VionicDentResin.

Storage/disposal

Storage at temperature, 15°C (59°F) to 28°C (82°F). Protect from direct sunlight. Even a small amount of light can trigger polymerization. To protect against contamination, cover the material in the tray with a lid or a glass plate. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Adverse effects

Product may cause allergic reactions.

Safety information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- When handling uncovered, use equipment with local exhaust ventilation.
- Wear suitable personal protective equipment (when handling the uncured material: protective gloves and goggles, when finishing the cured material additionally face mask).
- After contact with skin wash immediately with water and soap.
- If resin gets into your eyes immediately rinse and consult a doctor.
- Biocompatibility is only guaranteed with complete polymerisation.
- Please take note of the information in the latest version of the safety data sheet.
- Any serious incident that has occurred in relation to this product must be reported immediately to info@vita-zahnfabrik.com and to the competent authority of the Member State in which the user and/or patient is established.

MRI Safety Information: VITA VIONIC DENT RESIN is MR Safe as it is composed of materials that are electrically non conductive, nonmetallic and nonmagnetic.

Technical Data

Chemical composition	Wt%	Physical properties*	Value
Polyfunctional acrylates	45-80%	Flexural strength DIN EN ISO 10477*	> 100 MPa
Urethane bismethacrylate polymer	5-20 %	Flexural modulus DIN EN ISO 10477*	> 2800 MPa
Methacrylate monomer	0.1-5 %	Water sorption DIN EN ISO 10477*	< 40 µg/mm ³
Fillers (particle size 1.5 µm)	< 50%	Water solubility DIN EN ISO 10477**	< 7.5 µg/mm ³
Photoinitiators, pigments, stabilizers	< 0.2-10%	Viscosity**	1400-2200 mPas

*Applies to cured objects | according to norm at room temperature **Applies to liquid resin

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

Symbol Explanation

Graphic	Title	Description	Reference
	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements

Graphic	Title	Description	Reference
	Distributor	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Contains hazardous substances	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine-disrupting properties.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Medical device	Indicates that the product is a medical device.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Prescription only	Requires a prescription in the United States of America Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a physician.	USA Code of Federal Regulations 21CFR Part 801 § 801.109 (b)(1)
	Magnetic resonance	The freeprint crown is MR Safe.	FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021.
	Caution	The uncured material is considered as hazardous substance.	29 CFR 1910.1200 Occupational Safety and Health Standards
	This way up	This is the correct upright position of the distribution packages for transport and/or storage.	ISO 780:2015(E) Packaging-Distribution packaging- Graphical symbols for handling and storage of packages

VITA Zahnfabrik has been certified, and the product **VITA VIONIC® DENT RESIN** bears the mark:

CE 0483

R_{Only}

Made in Germany

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Spitalgasse 3

79713 Bad Säckingen

Germany

Phone: +49 7761 562-0

Fax: +49 7761 562-299

E-Mail: info@vita-zahnfabrik.com

detax GmbH

Carl-Zeiss-Str. 4

76275 Ettlingen

Germany

Phone: +49 7243 510-0

Fax: +49 7243 510-100

E-Mail: post@detax.com

922-13053_0125_V01

VITA